



Contact
medical innovations

EN ISO 13485:2018

MOBILE ENDOSURGERY SET ECONT M

The Mobile Endosurgical Set is designed for use in field hospitals operated by the Ministry of Defence or the State Emergency Service.

Highly mobile, the entire system is packed into two anti-shock water-proof transport containers, enabling endosurgical procedures and diagnostic examinations of the abdominal cavity that are usually only available in stationary clinical settings.

The Set includes:

- endoscopic video camera
- 24" TFT surgical video monitor
- endoscopic light source
- electrosurgical generator
- suction-irrigation unit
- electronic insufflator
- modular equipment trolley
- CO₂ gas supply module
- EcoFlow-type charging station
- fiber-optic light guides
- laparoscopes



Start-up time from transport to working condition - no more than 30 minutes.

Power supply for the set - 230 V, 50 Hz AC mains or EcoFlow unit.

Autonomous operation time from the EcoFlow unit - no less than 8 hours.

Weight in transport mode - no more than 120 kg.

Design and Manufacture:

Contact Co. LLC

104 Kyrylivska str.

04080 Kyiv, Ukraine

Ph. +380 44 490-9356

export@contact-endoscopy.com

<http://www.contact-endoscopy.com>